

Cost comparison of commercial SPRI bead options

Oct 2025


Omega and Illumina beads are ~50% the cost of Beckman beads.

Medium and Large sizes are much more economical than Small

Omega have been successfully integrated into library preparation workflows at UCSF (and [tested here](#)); Illumina IPB have been QC'd using 50bp ladder (see next slide), but not tested in the full workflow. We anticipate no issues.

Product	\$ per ml	Price \$USD	ml	Storage	Catalog number
Beckman AMPure XP Beads for DNA Cleanup	81.20	406	5	4C	A63880
Beckman SPRIselect DNA Size Selection Reagent	92.20	461	5	Room Temp	B23317
Omega Bio-Tek Mag-Bind® TotalPure NGS	25.20	126	5	4C	M1378-00
Illumina Purification Bead	N/A	Unavailable	10	Room Temp	20036089
Illumina Purification Bead	18.00	540	30	Room Temp	20119944
Beckman AMPure XP Beads for DNA Cleanup	25.07	1504	60	4C	A63881
Beckman SPRIselect DNA Size Selection Reagent	28.70	1722	60	Room Temp	B23318
Omega Bio-Tek Mag-Bind® TotalPure NGS	13.98	699	50	4C	M1378-01
Illumina Purification Bead	12.42	1242	100	Room Temp	20060057
Beckman AMPure XP Beads for DNA Cleanup	16.85	7581	450	4C	A63882
Beckman SPRIselect DNA Size Selection Reagent	18.60	8368	450	Room Temp	B23319
Omega Bio-Tek Mag-Bind® TotalPure NGS	8.54	4272	500	4C	M1378-02
Illumina Purification Bead	8.10	3240	400	Room Temp	20060058

Size Selection comparison between Illumina purification beads (IPB) and SPRISelect on a 50bp ladder



- IPB perform very similarly to SPRISelect, are also stored at room temperature, and cost much less (~50%)
- If there is a concern for specific sensitive applications, you may consider adjusting ratios up by 0.05% (for example if using 0.8x for SPRISelect, may consider 0.85x for IPB)
- We do not know how IPB perform after expiration. This has not been tested